

# Real-Time Soleus EMG Biofeedback for Gait Analysis and Stability Assessment

De Busschere Yévan – BIOMED

19/05/2025 - 25/07/2025

Irisse Laboratory - 117 rue du général Aileret - 97430 - Le Tampon - La Réunion

# Table of contents

1. Context and Motivation
2. Objectives of the Internship
3. Methods
  - a) Biofeedback System
  - b) Stability Assessment
4. Results
5. Discussion
6. Phelma Competence

# Context and motivation

Why do elderly people tend to walk more slowly, show reduced stability, and struggle to climb hills?

- Decrease in ankle push-off during walking
  - Increase in hip contribution

- Less efficient propulsion → slower speed, reduced stability
- Hip strategy → higher metabolic demand → faster fatigue

# Context and motivation

- **Common assumption** (Alfieri et al. 2010, Geirsdottir et al. 2012)  
Sarcopenia and leg muscle weakness → muscle atrophy → strengthening program  
→ do not improve the walking ability

- **Nonetheless**, Franz et al. (2013) proved that
  - Elderly adults (~72 yrs) exert 12% less propulsive force than young adults
  - With real-time visual feedback, they recovered or exceeded young adults' propulsion
  - **Reveals an unused ankle power reserve**



- Develop a real-time EMG biofeedback system
- Study its impact on stability

# Objective of the internship

1. Develop a real-time EMG biofeedback system targeting the soleus muscle
2. Provide visual feedback of muscle activation during push off phase
3. Assess whether modulating this activation affects whole body stability

# Background

- **Gait cycle** : focus on the terminal stance & pre-swing (push-off).
- **Soleus** : single-joint plantarflexor, main contributor to propulsion → ideal for feedback.
- **EMG** : measures muscle activation non-invasively
- **GRF** : Ground reaction forces measured by forces plates




Figure 1 :Calf anatomy




Figure 2: Raw EMG signal




Figure 4 : Gait cycle



Figure 3: Raw GRF (Ground reactions forces) signal

# Now that we have the objective and the background

1. How do we create a biofeedback?
  - a) Signal processing pipeline
  - b) System implementation (python)
  - c) Graphic User Interface
2. How is stability measured?
  - a) Motion capture system
  - b) Kinematic moment
3. Results

# How do we create a biofeedback ?



Figure 5 : Surface EMG




Figure 6 : Biofeedback setup



# How do we find the push-off phase?

- 50-Newton threshold in the vertical component → beginning of a stride
- Identify the positive phase of the longitudinal curve → compute mean EMG over that part




Figure 7 : Ground reactions forces

# Signal processing pipeline




Figure 8 : Signal processing pipeline

## Main Challenges :

- Achieve minimal delay in feedback
- Ensure real-time processing → *use lfilter instead of filtfilt*
- Maintain synchronization between signals

# System Implementation




Figure 9 : "Python pipeline"

- **Language:** Python (PyQt5 + qtm\_rt).
- **Multithreading** : separate threads for acquisition & GUI → ensure low latency.

# Graphic User Interface

- Main GUI : real time scatter of mean push-off EMG
- Feedback refreshes every stride
- Smoothed over 3 steps




Figure 10 : Main GUI

- Control window : Select condition (+/- 20%, 0%, define tolerance and save data)

The figure shows a 'Control Window' with various input fields and buttons. At the top, there are 'Start' and 'Enregistrer' buttons. Below them are two sliders: 'Consigne ±20%' and 'Tolérance ±10%', each with a blue arrow indicating the current setting. To the right of the sliders, there are buttons for 'Moyenne EMG', 'Gauche : 68.05', and 'Droite : 75.89'. Below the sliders, there are input fields for 'Masse du sujet (kg) : 80' and 'Seuil jambe gauche : 68', each with a 'Valider' button. At the bottom, there are three buttons labeled 'Condition 0%', 'Condition -20.0%', and 'Condition +20.0%'.

Figure 11 : Control window

# How is stability measured?

# How is stability measured?

Using whole-body kinetic moment

- Reflects how the body controls rotation around its center of mass.
- Stable gait = small, controlled angular momentum variations.



How do we measure the kinematic moment ?

# How is stability measured?

- Motion capture system ( marker + infrared cameras)
- Markers placed on the body to capture movement



Figure 12 : Treadmill of the lab



Figure 13 : Example of markers




Figure 14 : Marker set – positions of each markers ( Maldonado et al. )

# How do we assess the stability?




Figure 15 : Qualisys software window



Kinematic moment calculated using a Python program developed by PhD students

- Model reconstruction
- Extract data into a c3d file



Now that we've seen how the biofeedback was created and how stability is assessed, let's go through the experimental protocol and the main results.

# Experimental protocol




Figure 16 : Experimental protocol

# Results




Figure 17 : Kinetic moment along three axes with respect to normalized time

- Only 2 subjects only due to a lack of time → results not conclusive
- Slight increase in medio-lateral variability (Hy) in -20% condition → Possible link between reduced ankle activity and balance control.
- Kinetic moment calculated over the period where the condition was best respected (majority of points within the tolerance line)

# Discussion

- Real-time soleus feedback was technically feasible (stable in real time)
- Soleus = reliable feedback muscle for propulsion
- Preliminary data: modulating ankle activation may influence stability
- Limitations: small sample (n=2), treadmill vs. natural gait
- Future directions : EMG hip measurement, inertial units

# Competence Phelma




Figure 17 : Competence Radar – C1- Level 2

- Realization of the EMG biofeedback in real time
- Applied class knowledge for signal processing and Python implementation
- Used external resources to gain new skills in biomechanics, programming, and QTM
- Tested multiple system versions to achieve a functional setup
- Work documented through a detailed report and well-structured code